Product Instructions	English:	pages	2-14
Mode d'emploi	Français:	pages	15-30
Gebrauchsanweisung	Deutsch:	Seiten	31-45
Istruzioni per il prodotto	Italiano:	pagine	46-60
Instruções do Produto	Português:	páginas	61-75
Instrucciones	Español:	páginas	76-52

BinaxNOW[®]

Influenza A & B

Test Kit

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EMERGO EUROPE
P.O. BOX 18510
2502 EM THE HAGUE
THE NETHERLANDS

INTENDED USE

The BinaxNOW® Influenza A & B Test is an *in vitro* immunochromatographic assay for the qualitative detection of influenza A and B nucleoprotein antigens in nasopharyngeal (NP) swab and nasal wash/aspirate specimens. It is intended to aid in the rapid differential diagnosis of influenza A and B viral infections. Negative test results should be confirmed by cell culture.

Caution: Assay sensitivity was determined primarily using archived specimens. Users may wish to establish the sensitivity of this test on fresh samples.

SUMMARY AND EXPLANATION OF THE TEST

Influenza is a highly contagious, acute, viral infection of the respiratory tract. It is a communicable disease that is easily transmitted through the coughing and sneezing of aerosolized droplets containing live virus. Influenza outbreaks occur each year during the fall and winter months.¹ Type A viruses are typically more prevalent than type B viruses and are associated with most serious influenza epidemics, while Type B infections are usually more mild.

Rapid diagnosis of influenza A and B has become more important due to the availability of effective antiviral therapy. Rapid diagnosis of influenza can lead to reduced hospital stays, antimicrobial use and cost of hospital care.¹

The BinaxNOW® Influenza A & B Test provides a simple, rapid method for the diagnosis of influenza A and B using NP swab and nasal wash/aspirate specimens. The easy-to-use format and rapid results allow for its use in "STAT" testing where it can provide information to assist with treatment and hospitalization decisions.

PRINCIPLES OF THE PROCEDURE

The BinaxNOW® Influenza A & B Test is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to detect influenza type A and B nucleoprotein antigens in NP specimens. These antibodies and a control antibody are immobilized onto a membrane support as three distinct lines and combined with other reagents/pads to construct a test strip. This test strip is mounted inside a cardboard, book-shaped hinged test device.

Swab specimens require a sample preparation step, in which the sample is eluted off the swab into elution solution, saline or transport media. Nasal wash/aspirate samples require no preparation. Sample is added to the top of the test strip and the test device is closed. Test results are interpreted at 15 minutes based on the presence or absence of pink-to-purple colored Sample Lines. The blue Control Line turns pink in a valid assay.

REAGENTS AND MATERIALS

MATERIALS PROVIDED

Note: The materials provided in the test kit are sufficient for testing nasal wash/aspirate specimens only. If NP swab specimens will be tested, the Nasopharyngeal Swab Specimen Accessory Pack may be purchased.

BinaxNOW® INFLUENZA A & B TEST KIT

Test Devices: A cardboard, book-shaped hinged test device containing the test strip.

Transfer Pipettes: Fixed volume (100 µl) transfer pipettes used to transfer sample to the test devices. Use only pipettes provided by Binax or a calibrated pipette capable of delivering 100 µl sample volume.


Positive Control Swab: Inactivated influenza A and B dried onto swab.

Negative Control Swab: Inactivated *Streptococcus* Group A dried onto swab.


Elution Solution Vials for Control Swabs: Vials containing a fixed volume (0.5 ml) of elution solution used to prepare the Control Swabs for testing.


NASOPHARYNGEAL (NP) SWAB SPECIMEN ACCESSORY PACK

NP Swabs: Sterile foam swabs for use in the BinaxNOW® Influenza A & B Test. Other sterile flexible shaft NP swabs may be used in place of the Binax provided swabs. Refer to Specimen Collection and Handling section for details.

Elution Solution Vials for Swab Specimens: Vials containing a fixed volume (0.5 ml) of elution solution used to prepare the Control Swabs for testing. Transport media or saline may be used in place of Binax Elution Solution. Refer to Performance Characteristics – Transport Media section for details.

MATERIALS NOT PROVIDED

Clock, timer or stopwatch; nasal wash collection containers

PRECAUTIONS

1. For *in vitro* diagnostic use.
2. Leave test device sealed in its foil pouch until just before use.
3. Do not use kit past its expiration date.
4. Do not mix components from different kit lots.
5. The white sample pad at the top of the test strip contains reagents that extract the target antigen from the virus. To ensure optimum performance, add the sample **SLOWLY** (drop by drop) to the **MIDDLE** of this pad such that all of the sample volume absorbs into the pad.
6. Solutions used to make the control swabs are inactivated using standard methods. However, patient samples, controls, and test devices should be handled as though they could transmit disease. Observe established precautions against microbial hazards.
7. **INVALID RESULTS** can occur when an insufficient volume of specimen is added to the test device. To ensure delivery of an adequate volume, make certain that the lower shaft of the transfer pipette is full and does not contain air spaces before dispensing contents of the pipette onto the Sample Pad of the device. If air spaces are present, expel the specimen back into the container by squeezing the top bulb and redraw the specimen into the pipette. Use a new pipette if necessary.
8. When testing nasal wash/aspirate samples, avoid viscous areas of the sample when drawing specimen into the transfer pipette. If the pipette becomes clogged, such that the lower shaft of the pipette is not full, expel the specimen back into container by squeezing the top bulb and redraw the specimen into the pipette. Use a new pipette if necessary.
9. All transfer pipettes, elution ampoules and test vials are single use items – do not use with multiple specimens.

STORAGE AND STABILITY

Store kit at room temperature (59-86°F, 15-30°C). The BinaxNOW® Influenza A & B Test kit and reagents are stable until the expiration dates marked on their outer packaging and containers.

SPECIMEN COLLECTION AND HANDLING

Use fresh NP swabs and nasal washes/aspirates for optimal test performance.

Collect nasal washes in standard containers. Test as soon as possible. Washes can be held at 2-8°C for up to 24 hours prior to testing in the BinaxNOW® test.

Use sterile cotton, rayon, foam or polyester flexible-shaft NP swabs to collect nasopharyngeal sample. Calcium alginate swabs are not recommended for use in this test. Elute sample within one hour of collection. Test as soon as possible. Eluted swab samples can be held at 2-8°C for up to 24 hours prior to testing in the BinaxNOW® test. If needed, transport sample at 2-8°C in a leak-proof container.

Allow samples to warm to room temperature before testing in the BinaxNOW® test. Swirl gently to mix before testing.

QUALITY CONTROL

Daily Quality Control:

The BinaxNOW® Influenza A & B test has built-in procedural controls. For daily quality control, Binax suggests that you record these controls for each test run.

Procedural Controls:

- A. An untested device has a blue line at the "Control" position. If the test flows and the reagents work, this blue line will always turn pink in a tested device.
- B. The clearing of background color from the result window is a negative background control. The background color in the window should be light pink to white within 15 minutes. Background color should not hinder reading of the test.

External Positive and Negative Controls:

Good laboratory practice suggests the use of positive and negative controls to ensure that:

- test reagents are working; and
- the test is correctly performed.

BinaxNOW® test kits contain Positive and Negative Control Swabs. These swabs will monitor the entire assay. Test these swabs once with each new test kit opened. Other controls may be tested in order to conform with:

- local, state and/or federal regulations;
- accrediting groups, and/or;
- your lab's standard Quality Control procedures.

Refer to NCCLS EP12-A and 42 CFR 493.1202(c) for guidance on proper QC practices (U.S. customers only).

If the correct control results are not obtained, do not report patient results. Contact Binax Technical Service during normal business hours (EST).

- Phone: 1-800-323-3199 or 1-207-772-3988
- Fax: 1-207-761-2074

SAMPLE PREPARATION PROCEDURE

Nasal Wash/Aspirate:

Nasal wash/aspirate samples do not need preparation. Go to Test Procedure.

Nasopharyngeal Swabs:

Elute swab in 0.5 to 3.0 ml of saline or transport media by vigorously rotating the swab in the liquid. Refer to Performance Characteristics – Transport Media section for details. Go to Test Procedure. If eluting swab in the Binax Elution Solution, follow the Swab Elution procedure below.

Swab (Control & Patient) Elution using Binax Elution Solution:




1. The test kit contains test vials pre-filled with elution solution. Twist off the test vial cap.
2. Put the swab to be tested into test vial. Rotate the swab vigorously three (3) times in the liquid.
3. Press the swab against the side of the vial and turn as you remove it from the vial. This removes sample from the swab.
4. Discard the swab.
5. Test the liquid sample (from the test vial) in the BinaxNOW® test as soon as possible. Go to Test Procedure.

TEST PROCEDURE

1. Remove device from the pouch just prior to testing and lay flat on work bench.
2. Fill pipette by firmly squeezing the top bulb and placing pipette tip into sample. Release bulb while tip is still in sample. This will pull liquid into the pipette. Make sure there are no air spaces in the lower part of the pipette.
3. See arrow on test device to find White Sample Pad. **SLOWLY** (drop by drop) add entire contents of pipette (100 µl) to the **MIDDLE** of this pad by squeezing the top bulb.
4. Immediately peel off brown adhesive liner from the test device. Close and securely seal the device. Read result in window 15 minutes after closing the device. Results read before or after 15 minutes may be inaccurate.

Note: When reading test results, tilt the device to reduce glare on the result window, if necessary.

RESULT INTERPRETATION


For a **NEGATIVE SAMPLE**, the BLUE Control Line in the **BOTTOM THIRD** of the window turns a pink-to-purple color. No other line appears.

For a **FLU A POSITIVE SAMPLE**, the BLUE Control Line turns a pink-to-purple color AND a second pink-to-purple Sample Line appears above it in the **MIDDLE THIRD** of the window. Any Sample Line, even when very faint, is positive.

For a **FLU B POSITIVE SAMPLE**, the BLUE Control Line turns a pink-to-purple color AND a second pink-to-purple Sample Line appears above it in the **TOP THIRD** of the window. Any

Sample Line, even when very faint, is positive.

A test is **INVALID** if the Control Line remains BLUE or is not present at all, whether a Sample Line(s) is present or not. Repeat Invalid tests with a new test device. Call Binax Technical Service if the problem persists.

800-323-3199 or 207-772-3988

REPORTING OF RESULTS

Result	Suggested Report
Positive for Flu A	Positive for Flu A protein antigen.
Positive for Flu B	Positive for Flu B protein antigen.
Negative	Negative for Flu A and Flu B protein antigens. Infection due to Flu A and Flu B cannot be ruled out. Flu A and/or Flu B antigen in the sample may be below the detection limit of the test. Binax suggests culture of negative samples.

LIMITATIONS

A negative test result does not exclude infection with influenza A and B. Therefore, the results obtained with the BinaxNOW® Influenza A & B Test should be used in conjunction with clinical findings to make an accurate diagnosis.

The BinaxNOW® Influenza A & B Test detects both viable and non-viable influenza A and B. Test performance depends on antigen load in the specimen and may not correlate with cell culture performed on the same specimen.

Performance of the BinaxNOW® Influenza A & B Test has not been established for monitoring antiviral treatment of influenza.

Visibly bloody samples may not be appropriate for use in the BinaxNOW® Influenza A & B Test.

EXPECTED VALUES

The prevalence of influenza varies from year to year, with outbreaks typically occurring during the fall and winter months.¹ The rate of positivity found in influenza testing is dependent on many factors including the method of specimen collection, the test method used, geographic location, and the disease prevalence in specific localities. Type A viruses are typically associated with most serious influenza epidemics, while Type B are typically milder. In a multi-center clinical study conducted by Binax in the U.S. during the 2002 influen-

za season, the average prevalence of influenza A (as determined by viral cell culture) was 26% in nasal wash samples and 20% in NP swab samples. The average prevalence of influenza B was 21% in nasal wash samples and 20% in NP swab samples.

PERFORMANCE CHARACTERISTICS

The sensitivity and specificity of the BinaxNOW® Influenza A & B Test are equivalent to that of the individual Binax NOW® Flu A and Binax NOW® Flu B Tests as demonstrated by comparison studies using retrospective frozen clinical specimens and inactivated viral standards. The clinical performance of the individual Binax NOW® Flu A and Flu B tests versus conventional culture methods was originally established in a multi-center prospective study conducted during the 2002 flu season. The specificity of the BinaxNOW® Influenza A & B Test, as compared to cell culture/DFA, is equivalent to that of the individual BinaxNOW® Flu A and Binax NOW® Flu B Tests as demonstrated by a prospective study using fresh clinical specimens.

Clinical Studies:

BinaxNOW® Influenza A & B Test Performance vs. Cell Culture / DFA – Prospective Study

The performance of the BinaxNOW® Influenza A & B Test was compared to cell culture and/or DFA, and to the Binax NOW® Flu A Test and the Binax NOW® Flu B Test, in a prospective study conducted in 2004 outside the US. Nasopharyngeal (NP) swab and nasal wash / aspirate specimens, collected at multiple sites from children (less than 18 years of age) and adults (18 years or older) presenting with influenza-like symptoms, were evaluated in the Binax test at a central testing lab.

Forty-seven percent (47%) of the population tested was male, 53% female, 40% pediatric (< 18 years), and 60% adult (\geq 18 years). No differences in test performance were observed based on patient age or gender. There were no invalid tests reported.

One hundred and thirteen (113) NP swab specimens and 1 wash/aspirate specimen were tested. One hundred and eight (108) of the 114 samples tested were influenza negative by culture/DFA, and 6 samples were influenza positive. When compared to culture/DFA, the BinaxNOW® Test was 75% (3/4) sensitive and 100% (110/110) specific for detection of influenza A and 50% (1/2) sensitive and 100% (112/112) specific for detection of influenza B. There was 100% agreement between the BinaxNOW® Influenza A & B Test and the individual Flu A and Flu B Tests.

BinaxNOW® A & B Test specificity by sample type versus cell culture / DFA, including 95% confidence intervals, is listed below.

BinaxNOW® Influenza A & B Test Specificity vs. Cell Culture/DFA

Sample	FLU A SPECIFICITY				FLU B SPECIFICITY			
	-/-	+/-	% Spec	95% CI	-/-	+/-	% Spec	95% CI
NP Swab	109	0	100%	97-100%	111	0	100%	97-100%
Wash/Aspirate	1	0	100%	16-99%	1	0	100%	16-99%
Overall	110	0	100%	97-100%	112	0	100%	97-100%

BinaxNOW® Influenza A & B Test Performance vs. Binax NOW® Flu A and Flu B Tests:

Performance of the BinaxNOW® Influenza A & B Test was compared to the current NOW® Flu A Test on 306 retrospective frozen clinical samples and to the NOW® Flu B Test on 303 retrospective frozen clinical samples. All clinical samples were collected from symptomatic patients at multiple physician offices, clinics and hospitals located in the Southern, Northeastern and Midwestern regions of the United States and from one hospital in Sweden. Fifty-three percent (53%) of the population tested was male, 47% female, 64% pediatric (<18 years) and 36% adult (≥ 18 years). Nasal wash/aspirate specimens comprised approximately 57% of the samples tested, while NP swabs represented 42%. No differences in test performance were observed based on patient age and gender or based on sample type tested.

The BinaxNOW® Influenza A & B Test was 100% sensitive and 96% specific for detection of influenza A vs. the NOW® Flu A Test and 93% sensitive and 97% specific for detection of influenza B vs. the NOW® Flu B Test. Test performance by virus type (A vs. B), by sample type (swab vs. wash/aspirate), and overall, including 95% confidence intervals, is detailed in the following tables.

BinaxNOW® Influenza A & B Test vs. NOW® Flu A Test for Detection of Influenza A

Sample	FLU A SENSITIVITY				FLU A SPECIFICITY			
	+/+	-/+	% Sens	95% CI	-/-	+/-	% Spec	95% CI
NP Swab	30	0	100%	89-100%	96	1	99%	94-100%
Wash/Aspirate	47	0	100%	93-100%	123	9	93%	88-96%
Overall	77	0	100%	95-100%	219	10	96%	88-96%

BinaxNOW® Influenza A & B Test vs. NOW® Flu B Test for Detection of Influenza B

Sample	FLU B SENSITIVITY				FLU B SPECIFICITY			
	+/+	-/+	% Sens	95% CI	-/-	+/-	% Spec	95% CI
NP Swab	2	0	100%	29-99%	126	1	99%	96-100%
Wash/Aspirate	12	1	92%	66-98%	152	9	94%	90-97%
Overall	14	1	93%	70-98%	278	10	97%	94-98%

Binax NOW® Flu A and Flu B Test Performance vs. Cell Culture

Performance of the Binax NOW® Flu A and Flu B Tests was compared to cell culture on 373 prospective clinical samples collected as part of a multi-center study conducted during the 2002 Flu season at physician offices and clinics located in the Western and mid-Atlantic United States. Fifty-four percent (54%) of the population tested was male, 46% female, 90% pediatric (< 18 years) and 10% adult (≥ 18 years). Nasal wash/aspirates comprised 51% of the samples tested, while NP swabs represented 49%. No differences in performance were observed based on patient age and gender or based on sample type tested.

The Binax NOW® Flu A Test was 80% sensitive and 93% specific while the Binax NOW® Flu B Test was 65% sensitive and 97% specific when compared to cell culture. The performance of the two tests by sample type (swab vs. wash/aspirate) and overall, including 95% confidence intervals, is detailed in the following tables.

Binax NOW® Flu A Test vs. Cell Culture

Sample	FLU A SENSITIVITY				FLU A SPECIFICITY			
	+/+	-/+	% Sens	95% CI	-/-	+/-	% Spec	95% CI
NP Swab	29	8	78%	62-88%	133	12	92%	86-95%
Wash/Aspirate	40	9	82%	69-90%	133	9	94%	89-97%
Overall	69	17	80%	71-87%	266	21	93%	89-95%

Binax NOW® Flu B Test vs. Cell Culture

Sample	FLU B SENSITIVITY				FLU B SPECIFICITY			
	+/+	-/+	% Sens	95% CI	-/-	+/-	% Spec	95% CI
NP Swab	21	15	58%	42-73%	142	4	97%	93-99%
Wash/Aspirate	29	12	71%	56-83%	146	4	97%	93-99%
Overall	50	27	65%	54-75%	288	8	97%	95-99%

Analytical Sensitivity:

The BinaxNOW® test limit of detection (LOD), defined as the concentration of influenza virus that produces positive BinaxNOW® test results approximately 95% of the time, was identified by evaluating different concentrations of inactivated Flu A/Beijing and inactivated Flu B/Harbin in the BinaxNOW® test.

Twelve (12) different operators each interpreted 2 devices run at each concentration for a total of 24 determinations per level. The following results identify a concentration of 1.03×10^2 ng/ml as the LOD for Flu A/Beijing and 6.05×10^1 ng/ml for Flu B/Harbin.

Flu A/Beijing		
Concentration (ng/ml)	# Detected	% Detected
1.03 x 10 ² (LOD)	23/24	96
5.60 x 10 ¹ (Cutoff)	*	50
3.27 x 10 ¹ (High Neg)	4/24	17
True Negative	0/24	0

Flu B/Harbin		
Concentration (ng/ml)	# Detected	% Detected
6.05 x 10 ¹ (LOD)	23/24	96
2.42 x 10 ¹ (Cutoff)	11/24	46
1.51 x 10 ¹ (High Neg)	6/24	25
True Negative	0/24	0

*Linear regression was used to calculate a line equation, which was then used to project the cutoff concentration of Flu A/Beijing.

To demonstrate comparable analytical sensitivity of the BinaxNOW[®] Influenza A & B Test and the individual NOW[®] Flu A and Flu B Tests, the Flu A and B cutoff levels identified above were evaluated in the NOW[®] Flu A and Flu B Tests.

The A/Beijing cutoff sample detected 50% of the time in the BinaxNOW[®] Influenza A & B Test was also detected 50% (12/24) of the time in the NOW[®] Flu A Test when tested by six (6) operators interpreting a total of 24 devices. Likewise, the B/Harbin cutoff sample detected 46% of the time in the BinaxNOW[®] Influenza A & B Test was detected 10% (4/40) of the time in the NOW[®] Flu B Test when tested by ten (10) operators interpreting 40 devices.

These data demonstrate that the analytical sensitivity of the BinaxNOW[®] Influenza A & B Test is equivalent or better than that of the individual NOW[®] Flu A and B Tests.

Analytical Reactivity:

The seven (7) live influenza A strains and five (5) live influenza B strains listed below tested positive in the BinaxNOW[®] Influenza A & B Test at concentrations ranging from 10²-10⁶ CEID₅₀/ml. Although the specific influenza strains causing infection in humans can vary year to year, all contain the conserved nucleoproteins targeted by the BinaxNOW[®] test.²

Influenza Strain	ATCC #
Flu A/WS/33	VR-825
Flu A/NWS/33	VR-219
Flu A/Hong Kong/8/68	VR-544
Flu A/Aichi/2/68	VR-547
Flu A/New Jersey/8/76	VR-897
Flu A/Mal/302/54	VR-98
Flu A/Port Chalmers/1/73	VR-810
Flu B/Lee/40	VR-101
Flu B/Brigit	VR-786
Flu B/Russia/69	VR-790
Flu B/Hong Kong/5/72	VR-791
Flu B/R75	VR-789

Analytical Specificity (Cross-Reactivity):

To determine the analytical specificity of the BinaxNOW® Influenza A & B Test, 36 commensal and pathogenic microorganisms (27 bacteria, 8 viruses and 1 yeast) that may be present in the nasal cavity or nasopharynx were tested. All of the following microorganisms were negative when tested at concentrations ranging from 10^4 to 10^8 TCID₅₀/ml (viruses), 10^7 to 10^8 organisms/ml (bacteria) and 10^6 organisms/ml (yeast).

Bacteria	Viruses	Yeast
<i>Acinetobacter</i>	Adenovirus	<i>Candida albicans</i>
<i>Bordetella pertussis</i>	Coronavirus	
<i>Enterococcus faecalis</i>	Coxsackie B4	
<i>Escherichia coli</i>	Cytomegalovirus (CMV)	
<i>Gardnerella vaginalis</i>	Parainfluenza 1	
<i>Haemophilus influenzae</i>	Parainfluenza 2	
<i>Klebsiella pneumoniae</i>	Parainfluenza 3	
<i>Lactobacillus casei</i>	Respiratory Syncytial Virus (RSV)	
<i>Legionella pneumophila</i>		
<i>Listeria monocytogenes</i>		
<i>Moraxella catarrhalis</i>		
<i>Neisseria gonorrhoeae</i>		
<i>Neisseria meningitidis</i>		
<i>Neisseria sicca</i>		
<i>Neisseria subflava</i>		
<i>Proteus vulgaris</i>		
<i>Pseudomonas aeruginosa</i>		
<i>Serratia marcescens</i>		
<i>Staphylococcus aureus</i>		
<i>Staphylococcus aureus</i> (Cowan protein A producing strain)		
<i>Staphylococcus epidermidis</i>		
<i>Streptococcus</i> , Group A		
<i>Streptococcus</i> , Group B		
<i>Streptococcus</i> , Group C		
<i>Streptococcus</i> , Group F		
<i>Streptococcus mutans</i>		
<i>Streptococcus pneumoniae</i>		

Interfering Substances:

The following substances, naturally present in respiratory specimens or that may be artificially introduced into the nasal cavity or nasopharynx, were evaluated in the BinaxNOW® Influenza A & B Test at the concentrations listed and were found not to affect test performance. Whole blood (1%) did not interfere with the interpretation of negative BinaxNOW® test results, but did interfere with the interpretation of Flu A LOD positive samples. Therefore, visibly bloody samples may not be appropriate for use in this test.

Substance	Concentration
1 OTC mouthwash	20%
3 OTC nasal sprays	15%
3 OTC throat drops	15%
2 OTC throat sprays	20%
4-acetamidophenol	10 mg/ml
Acetylsalicylic acid	15 mg/ml
Albuterol	20 mg/ml
Chlorpheniramine	5 mg/ml
Dextromethorphan	10 mg/ml
Diphenhydramine	5 mg/ml
Guaiacol glycerol ether	20 mg/ml
Oxymetazoline	0.05%
Phenylephrine	50 mg/ml
Phenylpropanolamine	20 mg/ml
Rebetol	500 ng/ml
Relenza	20 mg/ml
Rimantadine	500 ng/ml
Synagis	0.1 mg/ml
Tamiflu	50 mg/ml

Transport Media:

The following transport media were tested in the BinaxNOW® Influenza A & B Test as negative samples (no virus present) and after inoculation with the LOD levels of Influenza A & B. Media did not impact BinaxNOW® test performance, with the media alone testing negative in the NOW® test and media inoculated with LOD Influenza A & B testing positive on the appropriate test line in BinaxNOW® test.

Amies Media
Hank's Balanced Salt Solution
M4 Media
M4-RT Media
M5 Media
Stuart's Media
Saline

Reproducibility Study:

A blind study of the BinaxNOW® Influenza A & B Test was conducted at 3 separate sites using panels of blind coded specimens containing negative, low positive, and moderate positive samples. Participants tested each sample multiple times on 3 different days. There was 97% (242/250) agreement with expected test results, with no significant differences within run (replicates tested by one operator), between run (3 different days), between sites (3 sites), or between operators (6 operators).

REFERENCES

- 1) Williams, KM, Jackson MA, Hamilton M. Rapid Diagnostic Testing for URIs in Children: Impact on Physician Decision Making and Cost. *Infect. Med.* 19(3): 109-111, 2002.
- 2) Dowdle, W.R, Kendal, A.P., and Noble, G.R. 1980. Influenza Virus, p 836-884. *Manual of Clinical Microbiology*, 3rd edition, In Lennette, et. Al (ed.). American Society for Microbiology, Washington, D.C.

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Souche grippale	ATCC #
Flu A/WS/33	VR-825
Flu A/NWS/33	VR-219
Flu A/Hong Kong/8/68	VR-544
Flu A/Aichi/2/68	VR-547
Flu A/New Jersey/8/76	VR-897
Flu A/Mal/302/54	VR-98
Flu A/Port Chalmers/1/73	VR-810
Flu B/Lee/40	VR-101
Flu B/Brigit	VR-786
Flu B/Russia/69	VR-790
Flu B/Hong Kong/5/72	VR-791
Flu B/R75	VR-789

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RÉFÉRENCES

- 1) Williams, KM, Jackson MA, Hamilton M. Rapid Diagnostic Testing for URIs in Children: Impact on Physician Decision Making and Cost. *Infect. Med.* 19(3): 109-111, 2002.
- 2) Dowdle, W.R, Kendal, A.P., and Noble, G.R. 1980. Influenza Virus, p 836-884, *Manual of Clinical Microbiology*, 3rd edition, In Lennette, et. Al (ed.). American Society for Microbiology, Washington, D.C.

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© 2004 Binax, Inc.
Binax, Inc.